

# ORANGE JUICE CLEANING INSTRUCTIONS

1. Turn OFF water on/off valve (if fitted).
2. Fill a clean food quality bucket with hand hot tolerable water.
3. Disconnect BiB Connector from Orange Juice BIB and connect to cleaning adaptor.  
(Can use BiB Connector cut from used box)
4. Immerse QCD / cleaning adaptor in bucket of hot water and pull through at dispense points to remove all traces of Orange concentrate.
5. Empty bucket and fill with recommended cleaning solution.
6. Repeat No. 4 with cleaning solution.
7. Leave cleaning solution in line for recommended time.
8. Empty bucket and fill with fresh cold water.
9. Repeat No. 4 with cold water ensuring all traces of cleaning solution are removed.
10. Disconnect cleaning adaptor and reconnect the QCD to the Orange Juice BIB.
11. Pull through all dispense points until orange concentrate pulls through at each.
12. Turn ON water on/off valve.
13. Taste product to check quality
14. If product is not satisfactory contact BevTec Services on 01403 248988

**RECOMMEND CLEANING ORANGE JUICE EVERY TWO WEEKS**